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Elena Jiltsova

Neurosurgery department

DAG HAMMARSKJOLDS VAG

75185 UPPSALA

Sweden

<References: 97222956-FA and 97049464-FA>

09 July 2025

Dear Elena Jiltsova,

Enclosed are two separate related notifications:

1- DBS IPG Feedthrough Wire Breaks (97222956-FA)

2- SureTek Burr Hole Cover Kit, part of the Vercise DBS System (97049464-FA)

Instructions:

1- Please read carefully the two (2) Field Safety Notices attached.

2- Complete and sign the enclosed Acknowledgement form.

3- Return the Form to Boston Scientific at **Email: NordicProdQuality@bsci.com** by **30 July 2025**.

Should you have any questions or require assistance regarding this matter, please contact your Boston Scientific representative.

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Attachments: - Acknowledgement Form
- Field Safety Notices

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09 July 2025

Urgent Field Safety Notice

Vercise Genus™ Deep Brain Stimulation (DBS) Implantable Pulse Generators (IPGs) with Feedthrough (FT) Wire Break(s) - IFU Update

Dear Elena Jiltsova,

In August 2024, Boston Scientific reminded DBS customers to follow the steps outlined in the labeling/IFU (Instructions for Use) to implant Vercise DBS IPGs within a subcutaneous pocket. Boston Scientific is now providing an update to the previous communication. An IFU revision is being implemented to include an additional warning statement.

“Warning: Submuscular implant in the pectoral region (i.e. subpectoral) may expose the Vercise Genus rechargeable IPG to frequent muscle tension forces which may cause component damage (feedthrough wire(s)) resulting in device malfunction that could lead to adverse events including explant of the IPG.”

Since launch of the Vercise Genus DBS System in 2020, Boston Scientific has received 11 complaints related to this issue.

The Vercise Genus DBS rechargeable IPG continues to meet safety and performance expectations when used in accordance with device labeling.

Our records indicate that your facility received some of the concerned product. **The table below provides a complete list of all affected products**, including Product Description and Material Number (UPN). Please note that **only the devices listed below are affected**. No other Boston Scientific product is involved in this Field Safety Notice and no product is being recalled.

Product Description	UPN
Vercise Genus™ R16 Implantable Pulse Generator Kit	M365DB12160
Vercise Genus™ R32 Implantable Pulse Generator Kit	M365DB12320

Surgical Implant Manual Document Number
92691261
92495783
92920693

Device Description

The Boston Scientific Deep Brain Stimulation (DBS) System provides a reversible therapy where structures in the brain are stimulated with small electrical pulses. A Boston Scientific DBS System utilizes current steering (also known as multiple independent current control or MICC) across eight or sixteen Contacts per DBS Lead to provide precise positioning of stimulation.

Recommendations and Instructions

1- Review the Surgical Implant Manual updates provided in **Appendix 1**. These updates will be found in the IFU once implemented.

2- If you are a facility that has sent products to another hospital or a facility within your network, ensure this notification is forwarded to them.

3- A patient letter is available to Healthcare Professional on demand which can be shared with the patient and/or included within the patient's medical record.

4- Monitor per relevant IFU recommendations for any clinical observations of high monopolar impedances, undesired sensation, sudden loss of therapy, return of pre-implant symptoms and/or Bluetooth connectivity challenges, as these may be signals of potential feedthrough wire break(s) (component damage).

5- To provide awareness of this information, share this letter with any other clinicians in your hospital who use the Boston Scientific DBS System.

6- Maintain a copy of this letter in your facility's records.

7- Continue to report all device-related incidents or quality concerns experienced with the use of these devices to Boston Scientific at BSN.ComplaintCallCenter@bsci.com.

8- **Please complete the attached Acknowledgement Form even if you do not have any affected product.**

9- **When completed, please return the Acknowledgement Form to your Boston Scientific office for the attention of Email: NordicProdQuality@bsci.com on or before 30 July 2025.**


Additional Information

Any adverse events or quality concerns associated with use of this product should be reported to Boston Scientific via email at BSN.ComplaintCallCenter@bsci.com

Your Competent Authority is being notified of this Field Safety Notice.

Patient safety is Boston Scientific's highest priority. We are committed to transparent communication with physicians and healthcare professionals to ensure you have timely, relevant product information for managing your patients and their devices. If you require additional assistance or more information regarding this communication, please contact your local Boston Scientific representative.

Yours sincerely,

John Donohue
Vice President, Quality Assurance
Boston Scientific

Appendix 1 – Vercise™ Deep Brain Stimulation Systems Surgical Implant Manual Updates

Table 1 below provides the planned updates to the Surgical Implant Manual, highlighted in blue pending review by relevant regulatory bodies.

Surgical Implant Manual Section	Planned Updates
Create the IPG Pocket and Tunnel the Lead Extension	1- Create a pocket for the IPG under the skin in a location that is either in the chest or in the abdomen on the same side of the patient as the DBS Lead(s) and Lead Extension(s) connection: a- Mark the location of the IPG pocket. b- Use the IPG template provided to outline the intended pocket to guide the optimal pocket sizing. Note: It is important to keep the pocket small to prevent the IPG from turning over. c- Non-rechargeable IPGs should be implanted no deeper than 2.5 cm. Rechargeable IPGs should be implanted no deeper than 2 cm. Communication, including device programming, could become ineffective at depths greater than 2.5 cm. IPG charging could become ineffective at depths shallower than 0.5 cm or greater than 2 cm (rechargeable IPG only). Note: The 2.5 cm depth restriction does not apply to Vercise PC IPG (DB-1140-S). The depth restriction does apply to all other non-rechargeable Boston Scientific IPGs. Note: For full body MRI scan eligibility, confirm that the IPG is implanted according to the instructions contained in the ImageReady™ MRI Guidelines for Boston Scientific DBS Systems. Warning: Submuscular implant in the pectoral region (i.e. subpectoral) may expose the Vercise Genus rechargeable IPG to frequent muscle tension forces which may cause component damage (feedthrough wire(s)) resulting in device malfunction that could lead to adverse events including explant of the IPG.

Surgical Implant Manual Section	Planned Updates
Implanting the IPG	15- Place the IPG in the subcutaneous pocket with the etched writing “This Side Up” facing the skin and parallel to the skin surface. Note: Make the pocket no deeper than 2 cm for rechargeable IPGs and no deeper than 2.5 cm for Vercise Genus non-rechargeable IPGs. Communication, including device programming, could become ineffective at depths greater than 2.5 cm. For rechargeable IPGs, IPG charging could become ineffective at depths shallower than 0.5 cm or greater than 2 cm. Warning: Submuscular implant in the pectoral region (i.e. subpectoral) may expose the Vercise Genus rechargeable IPG to frequent muscle tension forces which may cause component damage (feedthrough wire(s)) resulting in device malfunction that could lead to adverse events including explant of the IPG. Warning: Incorrect placement of the IPG in the pocket could require a revision surgery. a- Coil the excess Lead Extension length under or around the IPG perimeter. Warning: Avoid placing the excess length of the Lead Extensions on the superficial surface of the IPG as this may result in tissue erosion, ineffective communication, or charging difficulty. b- Optional: Secure the IPG to the fascia by suturing through the holes in the IPG Header. Caution: Do not excessively bend the Lead bodies and Connectors while placing the DBS System under the skin and closing the incision.

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09 July 2025

Urgent Field Safety Notice

SureTek™ Burr Hole Cover Kit, part of the Vercise™ Deep Brain Stimulation System, Retaining Clip Locking Difficulties – IFU Update

Dear Elena Jiltsova,

Boston Scientific is initiating a Product Correction (IFU update) to address difficulties while locking SureTek™ Burr Hole Covers (BHC) to secure Deep Brain Stimulation (DBS) Leads, part of the Vercise™ Deep Brain Stimulation System. In some cases, resistance has been encountered while closing the Retaining Clip locking mechanism to securely hold the lead in place.

This letter provides labelling updates to guide users if they encounter resistance or difficulties while closing the Retaining Clip locking mechanism. See **Appendix 1** for Surgical Implant Manual updates to be implemented.

The most common outcome is a limited procedural delay to exchange retaining clips or to withdraw the stylet before locking the Retaining Clip. If the lead is not properly secured during the initial procedure and there is resulting lead migration, the most serious outcome is the need for additional surgery to reposition a DBS lead.

Boston Scientific has received complaints describing prolonged procedures and/or the need for additional intervention due to difficulties using the SureTek Burr Hole Cover. Boston Scientific has received 46 complaints related to this issue, a rate of less than 0.2% over the past two-year period.

Our records indicate that your facility received some of the concerned product. **The table below provides a complete list of all affected products**, including Product Description, Material Number (UPN). **This letter only affects the surgical implant manuals listed below.** No other Material Numbers (UPNs) are impacted and no product is being recalled.

Product Name	UPN
SureTek™ Burr Hole Cover Kit	M365DB4600C0
	M365DB4605C0

Surgical Implant Manual Document Number
92691261
92495783
92920693

Product Description

The SureTek BHC consists of four implantable components referred to as the Base, Bone Screws, Retaining Clip, and Cap. The Base is positioned over a burr hole and secured with Bone Screws. The Retaining Clip is placed over the Base once a DBS lead is implanted. The Retaining Clip anchors the DBS lead with a locking mechanism enabled using the included Placement Tool.

The Surgical Implant Manual includes instructions for use of these devices.

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1- Review the updates to the Surgical Implant Manual provided in **Appendix 1**. These updates will be found in the IFU once implemented.

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3- To provide awareness of this information, share this letter with any other clinicians in your hospital who use the Boston Scientific Burr Hole Cover Kit or Burr Hole Cover Spares Kit.

4- If users encounter resistance or difficulties while closing the Retaining Clip locking mechanism, refer to the Surgical Implant Manual updates provided in **Appendix 1**.

5- Maintain a copy of this letter in your facility's records.

6- Continue to report all device-related incidents or quality concerns experienced with the use of these devices to Boston Scientific at BSN.ComplaintCallCenter@bsci.com

7- **Please complete the attached Acknowledgement Form even if you do not have any affected product.**

8- **When completed, please return the Acknowledgement Form to your Boston Scientific office for the attention of Email: NordicProdQuality@bsci.com on or before 30 July 2025.**

Additional Information

Any adverse events or quality concerns associated with use of this product should be reported to Boston Scientific via email at BSN.ComplaintCallCenter@bsci.com.

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Table 1: Updates to Instructions for Use

Surgical Implant Manual Section	Planned Updates
Securing the DBS Lead	6. While stabilizing the DBS Lead, carefully position the Retaining Clip over the Base so that the DBS Lead is located in the open channel of the Retaining Clip. Position the Retaining Clip so that the static side of the opening is against the Lead (Figure 21).  Figure 21. Position the Retaining Clip Over the Base 7. Push the Retaining Clip down into the Base. Ensure that the Retaining Clip is completely seated in the Base (Figure 22).  Figure 22. Push the Retaining Clip Into the Base 8. Place the tip end of the Placement/Removal Tool into the closure dimple or anywhere along the length of the Slider on the Retaining Clip to push the Slider towards the DBS Lead until it locks into place. Use the tip end of the Placement/Removal Tool to apply pressure on the Slider face in the opposite direction to ensure that the Slider is fully locked (Figure 23).

Note: If you have difficulty locking the Slider, you may need to remove the Lead Stylet before locking the Slider in place, or replace the Retaining Clip using another Burr Hole Cover Kit or Burr Hole Cover Spares Kit.


Figure 23. Lock the Slider

9. Remove the Lead Stylet. You may confirm the Slider is fully locked by repeating Step 8.

Caution: Do not reinsert the Lead Stylet into the DBS Lead while the DBS Lead is in the brain, as this may damage the DBS Lead and/or cause patient harm.

10. Verify that the DBS Lead has not moved from the desired location (e.g. by intraoperative imaging).

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